


# Blood Pressure Module

Accurate measurement without trigger

Original  
Quality  
Option



## Accurate measurement without trigger

An auscultatory blood pressure module with DKA (3-Dimensional K-sound Analysis) for use during ergometry stress testing. A complete pediatric cuff comes standard with this option. The microphone is double-sided.

### High standards

Lode is a socially and environmentally responsible company. All Lode products are RoHS/WEE compliant and Lode is ISO 9001:2015, and ISO 13485:2016 certified. All medical products comply to MDD 93/42/EEC, incl. IEC 60601-1.



### Original Accessory

Accessory designed and manufactured according to the Lode company quality standards



# Blood Pressure Module

Accurate measurement without trigger

Original  
Quality  
Option

## Specifications

### Accuracy

Minimum systolic measurement	50 mmHg
Maximum systolic measurement	250 mmHg
Minimum diastolic measurement	20 mmHg
Maximum diastolic measurement	150 mmHg

### Order info

Partnumber: 928820 -  
not  
available

*\*Specifications are subject to change without notice.*