

# Colour Display 3.5" - 2nd screen



Multifunctionality



## Multifunctionality

The display can be used on a treadmill or ergometer to display relevant data. Read out is available in 15 languages. For a treadmill: actual speed, actual elevation, testtime and Heart rate(when available) For an ergometer: actual RPM, actual power, time and Heart rate(when available) The colour display as a second screen can only be ordered in combination with a (programmable) 7 inch control unit. As a single screen, choose part number 945810).

### Multifunctional

The treadmill can be used in various ergometry fields, like cardiology, pulmonary function, rehabilitation, sports medicine and physiotherapy, as a stand alone unit or controlled by other equipment.

### Easy to operate

For Lode products this means:

- easy to connect
- easy to move around
- easy user interface

### High standards

Lode is a socially and environmentally responsible company. All Lode products are RoHS/WEE compliant and Lode is ISO 9001:2015, and ISO 13485:2016 certified. All medical products comply to MDD 93/42/EEC, incl. IEC 60601-1.



### Versatile Interfacing

Various interface protocols guarantee perfect communication with all commonly known stress ECG and spirometry equipment



### Robust design

The product is designed to withstand continuous heavy use by subjects in most weights and sizes.



### Original Accessory

Accessory designed and manufactured according to the Lode company quality standards



# Colour Display 3.5" - 2nd screen

Original  
Quality  
Option

Multifunctionality

## Specifications

### User Interface

English user interface

Chinese user interface

Croatian user interface

Czech user interface

Danish user interface

Dutch user interface

Finnish user interface

French user interface

German user interface

Greek user interface

Hungarian user interface

Italian user interface

Japanese user interface

Korean user interface

Latvian user interface

Lithuanian user interface

Norwegian user interface

Polish user interface

Portugese user interface

Romanian user interface

Russian user interface

Spanish user interface

Swedish user interface

Turkish user interface

Ukrainian user interface

Readout Angle

Readout Distance

Readout Speed

Readout Heartrate

Readout Time

Readout Mode

Screen size (diagonal)



8.9 cm

3.5 inch

### Order info

Partnumber: 945819

*\*Specifications are subject to change without notice.*



**FOR LIFE** | UNDERSTANDING  
MOVEMENT & PERFORMANCE

ISO 9001:2016 and ISO 13485:2016

Lode B.V.  
Zernikepark 16  
9747 AN Groningen  
The Netherlands  
Tel: +31 50 5712811  
Fax: +31 50 5716746  
E-mail: ask@lode.nl  
Internet: www.lode.nl