

# Colour display 3.5" - single display

Clear feedback

Original  
Quality  
Option



## Clear feedback

This screen allows for read out of the basic actual settings of the ergometer or treadmill. It is intended for use in a cardiopulmonary setup where the load or speed and elevation are set by the ECG or Pulmonary system and the ergometer or treadmill is used in terminal mode. Note that this display can not be used for operating the ergometer or treadmill. The language can be set to your liking.



# Colour display 3.5" - single display

Original  
Quality  
Option

## Clear feedback

The treadmill can be used in various ergometry fields, like cardiology, pulmonary function, rehabilitation, sports medicine and physiotherapy, as a stand alone unit or controlled by other equipment.

For Lode products this means:

- easy to connect
- easy to move around
- easy user interface

Lode is a socially and environmentally responsible company. All Lode products are RoHS/WEE compliant and Lode is ISO 9001:2015, and ISO 13485:2016 certified. All medical products comply to MDD 93/42/EEC, incl. IEC 60601-1.



## Versatile Interfacing

Various interface protocols guarantee perfect communication with all commonly known stress ECG and spirometry equipment



## Designed to be sweat-proof

The housing of the ergometer is designed in such way that sweat does not have the chance to drip into the mechanical parts. This ensures a long lifetime and makes the ergometer insensitive for malfunction.



## Original Accessory

Accessory designed and manufactured according to the Lode company quality standards

# Colour display 3.5" - single display

Clear feedback

Original  
Quality  
Option

## Specifications

### User Interface

English user interface	✓	
Chinese user interface	✓	
Croatian user interface	✓	
Czech user interface	✓	
Danish user interface	✓	
Dutch user interface	✓	
Finnish user interface	✓	
French user interface	✓	
German user interface	✓	
Greek user interface	✓	
Hungarian user interface	✓	
Italian user interface	✓	
Japanese user interface	✓	
Korean user interface	✓	
Latvian user interface	✓	
Lithuanian user interface	✓	
Norwegian user interface	✓	
Polish user interface	✓	
Portugese user interface	✓	
Romanian user interface	✓	
Russian user interface	✓	
Spanish user interface	✓	
Swedish user interface	✓	
Turkish user interface	✓	
Ukrainian user interface	✓	
Readout Angle	✓	
Readout Speed	✓	
Readout RPM	✓	
Readout Heartrate	✓	
Readout Time	✓	
Readout Power	✓	
Readout Mode	✓	
Screen size (diagonal)	8.9 cm	3.5 inch
Touchscreen	✓	

### Connectivity

Lode 38K4 interface protocol	✓
Lode ERM interface protocol	✓
Lode interface protocol	✓
Lode WLP interface protocol	✓
Ergoline P10 interface protocol	✓
Ergoline P4 interface protocol	✓
Schiller interface protocol	✓
"HP Cosmos" interface protocol	✓
"Woodway" interface protocol	✓
"Trackmaster" interface protocol	✓
"GE" interface protocol	✓

### Order info

Partnumber: 945810

*\*Specifications are subject to change without notice.*