

# Mobile safety belt and fall



## stop

Extra safety



### Highlights

#### High standards

Lode is a socially and environmentally responsible company. All Lode products are RoHS/WEE compliant and Lode is ISO 9001:2015, and ISO 13485:2016 certified. All medical products comply to MDD 93/42/EEC, incl. IEC 60601-1.

#### Safety for all users

The accessory enables safety lanyard use for subjects / patients with larger waist size

#### Connection with emergency stop

When the test subject falls, the treadmill will stop running

#### Mobile

The system can easily be moved

#### Comfortable Harness

The harness is very comfortable and does not obstruct the natural running movement

Harness that can be attached to a mobile safety device that prevents the test subject from falling on the belt. This construction on wheels can be placed over the treadmill whenever needed. The front wheels of the construction can be locked.

The height of the product is suitable for normal rooms (from 2,5 meters). For higher rooms you may consider to order height adjustment parts to enable use by taller patients. This will increase the height with 18 cm. Partnumber of this option is 932824.



# Mobile safety belt and fall stop

Extra safety

## Features



### Original Accessory

Accessory designed and manufactured according to the Lode company quality standards



### Comfortable harness

The comfortable harness is especially made for Lode to guarantee a comfortable fitting for the patient. The patient will not be burdened by the compensation weight, enabling the most comfortable way of rehabilitation. The harness is available in various sizes: small, standard and extra large.










# Mobile safety belt and fall stop



## Extra safety

Mobile safety belt and fall stop can a.o be extended with the following options:

<p><b>Harness Fallstop Small</b></p> <p>Suitable for smaller people</p>  <p>Partnumber: 458102</p>	<p><b>Harness Fallstop XL</b></p> <p>Suitable for larger people</p>  <p>Partnumber: 458105</p>	<p><b>Harness Fallstop XS</b></p> <p>Suitable for children</p>  <p>Partnumber: 458101</p>	<p><b>Harness Fallstop Large</b></p> <p>Suitable for larger people</p>  <p>Partnumber: 458104</p>	<p><b>Harness Fallstop Medium</b></p> <p>Standard supplied with Fallstop</p>  <p>Partnumber: 458103</p>
<p><b>Height adjustment 18 cm for BWSS / Mobile Fallstop</b></p> <p>Flexibility for tall patients</p>  <p>picture not yet available</p> <p>Partnumber: 932824</p>	<p><b>Height adjustment 36 cm for BWSS / Mobile Fallstop</b></p> <p>Flexibility for tall patients</p>  <p>picture not yet available</p> <p>Partnumber: 932833</p>			

# Mobile safety belt and fall stop



Extra safety

## Specifications

### Comfort

Allowed user weight 160 kg 352.7 lbs

### User Interface

English user interface	✓
Chinese user interface	✓
Croatian user interface	✓
Czech user interface	✓
Danish user interface	✓
Dutch user interface	✓
Finnish user interface	✓
French user interface	✓
German user interface	✓
Greek user interface	✓
Hungarian user interface	✓
Italian user interface	✓
Japanese user interface	✓
Korean user interface	✓
Latvian user interface	✓
Lithuanian user interface	✓
Norwegian user interface	✓
Polish user interface	✓
Portugese user interface	✓
Romanian user interface	✓
Russian user interface	✓
Spanish user interface	✓
Swedish user interface	✓
Turkish user interface	✓
Ukrainian user interface	✓

### Dimensions

Product length (cm)	160 cm	63 inch
Product width (cm)	150 cm	59.1 inch
Product height	250 cm	98.4 inch
Product weight	76 kg	167.6 lbs
Inner width	125 cm	49.2 inch

### Order info

Partnumber: 932845

*\*Specifications are subject to change without notice.*