

HUR by Lode



Evidence-based and smart rehabilitation

CARDIO with HUR SmartTouch ergometers and treadmills



Train and track the people you care for – your customers!

HUR SmartTouch is the perfect combination of training equipment and operating software, and is especially developed for senior exercise, rehabilitation and the inclusive wellness markets. Lode has developed a special range of ergometry products that is perfectly compatible with the HUR SmartTouch system. With this combination of quality products, you will enhance the training experience, you can operate wherever you need and whenever you want, you will be able to operate more efficient and save staff resources and you will simplify the work for the instructors.

CARDIO with HUR SmartTouch Treadmill

Modern designed treadmill specifically designed for rehabilitation. The low step-up height makes the treadmill perfectly suitable for all test subjects in the field of rehabilitation. The emergency stop with magnetic lanyard on the standard front handrail provides additional safety for the user. The treadmill offers a smooth acceleration from 0 km/h and is continuously adjustable in a range of 0,1 – 20 km/h.

Thanks to the unique low design of the cover plate it is possible to place a mirror or camera in front of the treadmill which makes it easier to monitor the gait of the test subject. The treadmill has a standard running surface of 50 x 150 cm, 0-25% elevation and heartrate. It is equipped with a 7" display with Touchscreen, a bottle holder and fixed handrails.



A USB A-B cable only for service purposes will be standard delivered with the product. To connect to the HUR smart touch server you only need a UTP network cable. The treadmill is standard executed with extended handrails.








HUR by Lode

Evidence-based and smart rehabilitation



Lode HUR by Lode product range consists of the following products:

<p>CARDIO with HUR SmartTouch Treadmill</p> <p>Robust and reliable treadmill for rehabilitation in the SmartTouch ...</p>  <p>Reliable treadmill with modern design</p>	<p>CARDIO with HUR SmartTouch Ergometer</p> <p>Optimal adjustability and low effort ergonomic bikemounting</p>  <p>Comfortable and reliable ergometer</p>	<p>CARDIO with HUR SmartTouch Recumbent Ergometer</p> <p>Easy positioning and entrance for all test subjects</p>  <p>Ergometer for elderly and obese test subject</p>	<p>CARDIO with HUR SmartTouch Arm Ergometer</p> <p>Modern ergometer with multifunctional applications</p>  <p>Smarttouch arm ergometer with automatic stand</p>	<p>CARDIO with HUR SmartTouch Arm Ergometer with Multifunctional Application</p> <p>Modern ergometer with multifunctional applications</p>  <p>Arm ergometer with adjustable wall fixation</p>
---	---	---	---	--

About Lode

Ever since Freerk Lode manufactured the first electro-magnetic bicycle ergometer in 1952, accuracy, reliability and durability have been fundamental for further developments. Having started years ago in the market of cardiology and pulmonary function, Lode BV has become a specialist in the complete spectrum of medical ergometry. Lode is world renowned as a manufacturer of high quality ergometers and the Lode brand stands for accuracy, durability and ergonomic design. The Lode product range varies from bicycle ergometers and treadmills to recumbent, arm and supine ergometers and ergometry software. Long-term experience in manufacturing medical equipment and continuous development to meet the changing requirements of the market, make Lode a flexible and reliable partner. Together we can transform your specific ideas and wishes into custom-made products. Before leaving the factory all Lode ergometers are dynamically calibrated and, of course, all units are produced under the strictest quality control conditions. Lode is ISO 9001:2015, ISO 13485:2016 certified and fulfils the EU Medical Device Directive MDD 93/42/EEC. Many Lode products are FDA registered. Over years of use, service costs are almost negligible.



ISO 9001:2015 and ISO 13485:2016

Lode B.V.
Zernikepark 16
9747 AN Groningen
The Netherlands
Tel: +31 50 5712811
Fax: +31 50 5716746
E-mail: ask@lode.nl
Internet: www.lode.nl